

Anesthesia Tech Step-by-Step: Awake Fiberoptic Intubation

- Confirm
 - **Weight** ____ kg, **BMI** ____ , and Dexmedetomidine ordered
- Retrieve...
 - Awake Fiberoptic Kit
 - 2% lidocaine x 2
 - Alaris Pump – 2 channel
 - Fiberoptic scope and tower
- Setup...
 - Fiberoptic tower on patient's **LEFT** side head of bed
 - 10ml attached to atomizer
 - (**Optional**- check with provider) 5ml attached to epidural catheter through suction port
 - Load lubricated Parker 7.0 ET tube
- Retrieve
 - Dexmedetomidine (5pm-6am go to 6th floor OR front desk)
- Program...
 - Alaris pump for Dexmedetomidine
 - BOLUS: Dexmedetomidine **1mcg/kg over 10 minutes** (Actual body weight), (requires faux infusion 0.4mcg/kg/min to select BOLUS option)
- Confirm...
 - Dexmedetomidine programming with provider
 - Plan for awake fiberoptic intubation with OR RN

PLEASE PROVIDE FEEDBACK

Was the protocol followed?

Did you require any additional supplies?

Additional comments:

Please also collect provider feedback and return to level 6 induction room AFOI feedback return bin. Thank you!